

# SIUI CTS-5500 & CTS-7700



- Portable and Compact
- Large 10 inch Monitor
- Programmable Application Settings
- Multi-Frequency Transducers
  - Calculation Packages
  - Hard Drive and USB
  - Annotation Library
- Two probe connectors
- 256 or 512 frame Cine
  - Economical

## Specifications

Scanning Method	Linear, Convex, Endocavity
Display Modes	B, Dual B, Quad B, B/M, M
Gray Scale	256
Focal Zones	up to 4 variable
Image Adjustments	Gain, TGC 8 slides, Dynamic Range, Edge Enhancement, Post Processing , Persistence
# Frequencies per Probe	5
Measurements	Distance, Area, Circumference, Volume, OB/GYN, Urology, Abd, Cardiac, Small Parts
Monitor	10 inch
Probe Connections	2
Dimensions	12.2 in.W , 12.6 in.H , 16.4 in.D .
Weight	Approx. 26.5 lbs.

## Probes and Applications

Convex 60	2.5/3.0/3.5/ 4.0/5.0 MHz.	Abd, OB/GYN
Convex 40	2.5/3.0/3.5/ 4.0/5.0 MHz.	Card, Abd, OB
MicroConvex	4.5/5.2/6.0/ 6.7/7.5 MHz.	Ped Card, Ped Abd
Linear 50mm 5500 only	5.0/6.0/7.5 8.5/10 MHz.	Mammography, Thyroid, MSK, Superficial
Linear 38mm	5.0/6.2/7.5 8.2/10 MHz.	Mammography, Thyroid, MSK, Superficial
Linear 25mm 7700 only	8/9/10/11/12 MHz.	Small Parts, MSK
Endocavity	3.5/4.0/5.0/ 6.0/7.5 MHz.	Transvaginal Transrectal

Options:  
Video Printer, Footswitch, 3D on 7700  
Biopsy Guides, Dicom, Cart



## DIAGNOSTIC HEALTH GROUP

PO Box 747 • Frederick, MD 21703 • E-Mail: [DHGsonomed@aol.com](mailto:DHGsonomed@aol.com)  
800 669-3442 • Fax 301 695-7940 • [www.dhgsonomed.com](http://www.dhgsonomed.com)