

The OB1 Obstetrical Doppler



- Large display of the Fetal Heart Rate updates continuously and is easy to read.
- Detachable waterproof transducer. Recessed transducer and cable holder for storage.
- Rechargeable NiCad batteries.
- Large speaker with adjustable volume provides good sound quality.
- Five year warranty.

Hokanson[®]

OB1

OB1 Specifications



Features:

Liquid Crystal Display (LCD) is large and easy to read. Displays and continuously updates the Fetal Heart Rate (FHR), while a symbol indicates each detected heart beat. Also has a low battery indicator, and volume level display.

Four buttons provide simple user controls. Waterproof transducer is detachable and can be gas sterilized at temperatures below 140°F.

Doppler Frequency: 2.25 MHz

Instrument Size:

Length: 20.3 cm (8 in)
Width: 5.0 cm (2 in)
Height: 12.7 cm (5 in)

Transducer Size:

Length: 10.4 cm (4.1 in)
Width: 2.8 cm (1.1 in)

Weight: 740 gm (1.63 lbs)

Power: Rechargeable NiCad batteries

Audio: 7.5 cm (3 in) loudspeaker
Optional headphones.

Warranty:

The OB1 includes the exceptional Hokanson five year warranty. Hokanson instruments are warranted against defects in parts and workmanship for a period of five years from the date of delivery. Transducers and accessories are warranted for one year.

Hokanson[®]

D.E. Hokanson, Inc.
12840 NE 21st Place
Bellevue, WA 98005 USA
(425) 882-1689 or
(800) 999-8251
fax (425) 881-1636
www.hokanson.cc