


Hokanson NIVP3 Software



- Windows® Based
- Color Reports
- Patient Database
- Complete Range of Testing
 - Several Versions
- Can be used with several Hokanson products

The NIVP3 software package allows you to use several Hokanson Vascular systems for arterial and venous exams. The software and computer become an extension of your diagnostic instrument enhancing the power of both. Patient reports are comprehensive, complete with Doppler and plethysmographic waveforms, patient demographics and tabular data. The software is Windows® based that captures and prints your studies. The software can be used to track a vascular patients changes over time.

Optional versions include:
Spectral Analysis
Arterial Inflow



DIAGNOSTIC HEALTH GROUP

PO Box 747 • Frederick, MD 21703 • E-Mail: DHGsonomed@aol.com
800 669-3442 Fax 301 695-7940