



# SHIMADZU

## SDU-350XL

- Compact
- Portable
- Programable User Settings
- Dual Frequency Transducers
  - Calculation Packages
    - PC Card Archival
  - Annotation Library
    - Easy to Use
  - Economical



### Specifications

Scanning Method	Linear, Convex
Display Modes	B, Dual B, B/M, M
Gray Scale	256
Focal Zones	1 or 2
Image Adjustments	Gain 0-40, STC/TGC 6 slides, Dynamic Range 45-80dB, Edge Enhancement 4 steps, Post Processing 8 steps, Persistence 6 steps
Presets	9
Measurements	Distance, Elipse, Area, Volume, OB package, GYN package, Urology package, Cardiac package, Thyroid
Monitor	9 inch
Probe Connections	1, 2 Optional
Dimensions	W 300mm 12in, H 275mm 11in, D 425mm 17in.
Weight	15 Kg. 32 lbs.

### Probes and Applications

Convex 57R	2.5-6 MHz.	Abd, OB
Convex 15R	2.5-6 MHz.	Card, Abd, OB
Convex 15R	3-7 MHz.	Ped Card, Ped Abd
Linear 60mm	5-9 MHz.	Mammography, Thyroid Superficial
Linear 34mm	6-12 MHz.	Mammography, Thyroid Superficial
TV11R	4-8 MHz.	Transvaginal
TR10R Bi-plane	4-8 MHz	Transrectal need Dual conn

Options:  
Video Printer, VCR, Footswitch,  
Dual Probe Connector, Carrying Case



## DIAGNOSTIC HEALTH GROUP

PO Box 747 • Frederick, MD 21703 • E-Mail: DHGsonomed@aol.com

800 669-3442 Fax 301 695-7940