

# Nicolet Versalab APM



- Portable
- Numeric or Graphic Display
  - Integrated Printer
  - Selectable Alarm Settings
- Mother & Clinician Event Markers
  - Battery or Line Operated
  - Economical

The Versalab APM is a portable, economical antepartum monitor (APM) for performing non-stress tests. Its small size and light weight make it the ultimate in portability and convenience. The integrated printer provides immediate, real time print-outs of the exam. No waiting for results or interpretation. This documentation is also necessary for reimbursement and can make placing the document in the correct files easier since it is available sooner. Maternal and Clinician event markers and large display are features that make this system easy to use. Since the system can be run on batteries or plugged into a wall outlet it is practical for office, clinic, portable or in-home use.

## Specifications

Weight - Main Unit	4.1 Lbs (1.9 Kg)
Dimensions	W 12.2 in. (31cm) D 9.1 in. (23cm) H 2.9 in. (7.3cm)
Doppler Technology	Continuous Wave (CW) unfocused 2 MHz probe
Battery Type	7.2V NiMH 6x4/5A rechargeable



---

---

## DIAGNOSTIC HEALTH GROUP

---

---

PO Box 747 • Frederick, MD 21703 • E-Mail: [DHGsonomed@aol.com](mailto:DHGsonomed@aol.com)  
800 669-3442 Fax 301 695-7940